

SUMO3-SMT3 suppressor of mif two 3 homolog 3 (*S. cerevisiae*) Gene

Reference: PTXBC000036

| | | | |
|----------------|---|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | SMT3A; SMT3H1; SUMO-3; Smt3B; small ubiquitin-related modifier 3; SMT3 suppressor of mif two 3 homolog 1; SMT3 suppressor of mif two 3 homolog 3; ubiquitin-like protein SMT3B; small ubiquitin-like modifier 3 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | SMT3 suppressor of mif two 3 homolog 3 (<i>S. cerevisiae</i>) |
| GENE ID | 6612 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC000036 |
| VECTOR | pDONR233 |
| NCBI SYMBOL | SUMO3 |
| LENGTH | 312 bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atgtccgaggagaagcccaaggagggtgtgaagacagagaatgaccacatcaacctgaagggtgcccgggcaggacggctcctgtgtgcagttcaagatcaaggacacacgccgctgagcaagctgatgaaggcctactcgagaggcaggg
ctgtcaatgaggcagatcagattcaggttcgacggcgagccaatcaatgaaactgacactccagcagctggagatggaggacaggacacatcgacgtgtccagcagacgggagggtgtccggagagcagcctggcagggcacagt
ttctag