

SUMO3-SMT3 suppressor of mif two 3 homolog 3 (S. cerevisiae) Gene

Reference: PTXBC000036

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMT3A; SMT3H1; SUMO-3; Smt3B; small ubiquitin-related modifier 3; SMT3 suppressor of mif two 3 homolog 1; SMT3 suppressor of mif two 3 homolog 3; ubiquitin-like protein SMT3B; small ubiquitin-like modifier 3		

For research use only.

GENE DESCRIPTION	SMT3 suppressor of mif two 3 homolog 3 (S. cerevisiae)
GENE ID	6612
ORIGIN SPECIES	Human
ACCESSIONS	BC000036
VECTOR	pDONR223
NCBI SYMBOL	SUMO3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgaggagaagcccaaggagggtgtgaagacagagaatgaccacatcaactgaaggtggccggcaggacggctccgtggtgcagttcaagatcaagaggcacacgccgctgagcaagctgatgaaggcctactcgagagggcagggcctgtcaatgaggcagatcagatcaggttcacggcgagccaatcaatgaaactgacctccagcacagctggagatggaggacgagacacatcgactgttccagcagcagcggagggtgtgccggagagcagcctggcagggcacagtcttag