

SPIN2B-spindlin family, member 2B Gene

Reference: PTXBC000044

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPIN-2; SPIN-2B; SPIN2_duplicate; TDRD26; dJ323P24.2; spindlin-2B; spindlin family, member 2 duplicate; spindlin-like protein 2B; telomeric SPIN2 copy; spindlin family member 2B		

For research use only.

GENE DESCRIPTION	spindlin family, member 2B
GENE ID	474343
ORIGIN SPECIES	Human
ACCESSIONS	BC000044
VECTOR	pDONR233
NCBI SYMBOL	SPIN2B
LENGTH	777 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgaagacccccaacgcacaggaagccgaaggcaacaaccaggcagctgcaggacggccactgggtctgcaaacatgacaaagaaaaagtctccaaaagaagcagagagccgacctcctccagccccgcaggaacatcgtgg
gctgcagaatttctcatggatggaagggaaggatgagccatcacgcagtggaaggaaaccgttctgcatcaggcgcctataaatccctctcttctggtgaaatgatggaaatgactgtgtctatggactggaaactcacagagatgaagggtt
tgtctctaaattctctgacaggggtgcatcatctcacattagtgatgccaacctgcaaataccataattgcaaaagcagtggaacatattttgaggagagcatggttctaaggatgaatggaggggatggcttagctcaagcacctatcatgaa
agcctggtttatattacctatgagaagatcctgtctgtacatgtaccagctctagatgattataaggaaggtagacctcccatcagcagaatccagtgagtcctccaacagagaggagccagggaggtgtagatgccctaataaggaagca
tgtggaatataccaaaagatgctcctcaaaaggatcgcatggtcattcaccaggtaggaaccaaacctctgtgtatttcaagaattgatgatgatttcataatctatgtctacgatttggtaaaaagtctaa
```