

ICAM4-intercellular adhesion molecule 4 (Landsteiner-Wiener blood group) Gene

Reference: PTXBC000046

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD242; intercellular adhesion molecule 4; CD242 antigen; LW antigen a; Landsteiner-Wiener blood group antigen a; Landsteiner-Wiener blood group glycoprotein; intercellular adhesion molecule 4 (LW blood group); intercellular adhesion molecule 4 (Landsteine		

For research use only.

GENE DESCRIPTION	intercellular adhesion molecule 4 (Landsteiner-Wiener blood group)
GENE ID	3386
ORIGIN SPECIES	Human
ACCESSIONS	BC000046
VECTOR	pDONR223
NCBI SYMBOL	ICAM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgggtctctgtccctctgtcgtcgtgttttttggcggccctaccgggagttgggagcgcgtgggacgccgactaagcgggcgcaagcccaagggtagccctctcgcgccctccggacctcagtcctctctgggtgcgatgcccggagttcggctgtgcag
ccggggaagtcagtcagctcaatfgcagcaacagctgtcccagccgcagaattccagcctccgaccccgctcggcgaaggcaagacgtcagaggccgggttgggtgtcttaccagctgtcgtcagctgagggcctggagctccctcgcgcactgcctcgtgacctgcgcagga
aaaacacgtcggccacctccagatcaccgctacagttcccggtgggctactgggtgacctgagcctgggaagccgggtcattatccgaaagcctggagcgttccaccggcctggatctggccaactgacctgacctacgagtttctgctggaccccgactctgg
cagccctgatctgccacgcgcctcaatctgacggcctgggtccgcaacagctcggcaccattacatgatgctcgttggagccccgcgccacagcttggcctccggttccatcgtgccttggatggatcctcctcactgtggcgctcgtacatgcaagtgcctagc
tatgaagtcaccagcgtaaaagggggatgttctatgcccgtgagcgagaaaaaggaatatgaacaatctgggaaatggccatactggtgctga
```