

DNAJC17-DnaJ (Hsp40) homolog, subfamily C, member 17 Gene

Reference: PTXBC000048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dnaJ homolog subfamily C member 17; DnaJ (Hsp40) homolog, subfamily C, member 17; DnaJ heat shock protein family (Hsp40) member C17		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 17
GENE ID	55192
ORIGIN SPECIES	Human
ACCESSIONS	BC000048
VECTOR	pDONR233
NCBI SYMBOL	DNAJC17
LENGTH	915 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggcagtgaccaaggagctcttacagatggacctgtacgcgctgctaggcattgaggagaaggcagcggacaagaaggtaagaaggcgtataggcagaaggccctctcctgccaccagacaaaaatccagataatcccagagcagctgaact
ctccaccagctttcagccctggagggtgctgaccgatgctgcagccagggtgcatatgacaaggctcaggaaagccaagaagcaagcagcagagaggaccagaacctgatgagaaaaaggaagaagtgaagctgacctggaggccgg
gagcggcagggcccaggcccaggagagtgaggaggaaagagagagccggagcaccaggactagagcaagagatgcaacgcctgagagaaggggttcccggcagctggaggaacagcagaggctcatccgggagcagatacggcagg
agcgtgaccagagggtgagaggaaaggcagaaaatactgaaggccaaggaacccccaaactaaagctaaaatggaagtcaagaaggagatgagtcaaaagggtgctactccaaagacgtcctcctacggctttgcagaagtatggtgaggtt
ctcaacctggtgctttccagtaagaagccaggcactgctgtgtggtgagttgcaaccgtcaaggcagcggagctggtgtccagaatgaagtggcctggtgataaccctctgaaatttctggtgagggacagccccagatgccgtgggccc
cagccactcaggactgtcaaaaggctcagctgtcagagaggactacgagaccctgtcatgatgctgcatgcccaggcggccgagcggcaacagctgatgcacggatgcagcaggaagaccaggaggggccgcctacgtag
```