

**DNAJC17-DnaJ (Hsp40) homolog, subfamily C, member 17 Gene**

Reference: PTXBC000048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dnaJ homolog subfamily C member 17; DnaJ (Hsp40) homolog, subfamily C, member 17; DnaJ heat shock protein family (Hsp40) member C17		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	DnaJ (Hsp40) homolog, subfamily C, member 17
<b>GENE ID</b>	55192
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC000048
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	DNAJC17
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

```
atggcagtgaccaaggagctctacagatggacctgtacgcgctgctagcattgaggagaagcagcggacaagaggtaagaaggcgtatagcagaaggccctctccaccagacaaaaatccagataatccagagcagctgaactttccaccagctttcagccctt
ggaggtgctgaccgatgctgcagccagggctgcatatgacaaggtcaggaaaagccaagaagcaagcagcagagagagaccagaaaactgatgagaaaaagaaagtgaaagctgacctggagcccgagcggcagccagccagagagtgaggagaaagga
gagccggagcaccaggacactagagcaagagatcgaacgctgagagaagaggggtcccgagctggaggaacagcagagggctcaccggagcagatacggcaggagctgaccagaggttgagagaaagcagaaaatactgaaggccaaggaaccccaactaaag
ctaaaatggaagtgcaagaaggagatgagtaaaaggtggctactccaagacgtctctacggctttgacagaagtatggtgaggtctcaactgggtcttccagtaagaagccaggcactgctgggtggagttgcaaccgtcaaggcagcggagctggcttccagaatgaag
ttggcctggtggataaccctggaagatttctgtggaggacagccccaggatccgtgggcccagccactcaggactgtcaaaaggctcagtgctgtcagagaggactacgagagcctgcatgatgctgcatgcccagggccgagcggccaacagctgatgcacgga
tgcagcaggaagaccaggagggccgctacgtag
```