

TMEM38B-transmembrane protein 38B Gene

Reference: PTXBC000049

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf87; D4Ert89e; OI14; TRIC-B; TRICB; bA219P18.1; trimeric intracellular cation channel type B; transmembrane protein 38B		

For research use only.

GENE DESCRIPTION	transmembrane protein 38B
GENE ID	55151
ORIGIN SPECIES	Human
ACCESSIONS	BC000049
VECTOR	pDONR233
NCBI SYMBOL	TMEM38B
LENGTH	876 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattcctccatgggacgagttggctctggccttccccgcagtcctatgtttccctttttgacatcgcgcactatctagtgtcagtgatggcggtgaaacgtcagccgggagcagctgcattggcatggaagaatcctattcaagctggttactgctat
getccactgtttggtggaggaattttactctgtctactgcttgcagagcctccatggaagtctctgcaaacacactaacataactgcatcttcaatctggatattacattttttgcccgcacatgacatggttccagggctattcattctactctgttcaact
actggcttgggaatgaaggagtgaccagaactggaaaatagtaggtggagtcacacatgctaatagtctattacaaaaatggctggatagtcagctatggatggcccgaggtgcaggtgtaccattataacgaatttgagagggtgta
aaaggagattggaaccagaaggtgatgaatggctgaagatgtacacctgccaaggtaacctgctgggtcagttatcttcacattccagcacaccagcactctggcaatcaaaagcataatcttatgttctttataccatctttattgtggccaaa
agataacatgatgactacacagactctactatgacattgctcttttggagatacattgagttggatgctatttggctggcagcagccgttttcatcatgtgagaagaaaagtgaagcaaaagtcacctccaatggcgttgggtcattggcctcaagccg
gtagatgttgcctcagataatgttaaaaaaacataactaagaagaatgaataa