

SNRPG-small nuclear ribonucleoprotein polypeptide G Gene

Reference: PTXBC000070

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMG; Sm-G; small nuclear ribonucleoprotein G; sm protein G; snRNP-G; small nuclear ribonucleoprotein polypeptide G		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein polypeptide G
GENE ID	6637
ORIGIN SPECIES	Human
ACCESSIONS	BC000070
VECTOR	pDONR223
NCBI SYMBOL	SNRPG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcaaaagctcaccctccgagttgaaaaatttatggacaagaagtatcattgaaatfaaatgtgacagacatgtccaaggaattgcggggattgatcccttataaaccttgatagatgaatgtggagatggcgactagtgacaacagacaatattggaatggtgtaata
cgaggaaatagtatcatcatgttagaagccttgaacgagtataa