

SNRPG-small nuclear ribonucleoprotein polypeptide G Gene

Reference: PTXBC000070

| | | | |
|----------------|--|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | SMG; Sm-G; small nuclear ribonucleoprotein G; sm protein G; snRNP-G; small nuclear ribonucleoprotein polypeptide G | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | small nuclear ribonucleoprotein polypeptide G |
| GENE ID | 6637 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC000070 |
| VECTOR | pDONR233 |
| NCBI SYMBOL | SNRPG |
| LENGTH | 231 bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atgagcaaaagctcacctcccgaattgaaaaaattatggacaagaagtatcattgaaattaaatggtggcagacatgtccaaggaatattcggggattgatccctttatgaaccttgtgatagatgaatgtgtggagatggcgactagtggacaacag
aacaatattggaatggtgtaatacagagaaatagatcatcatgttagaagccittgaacgagtataa