

## MRPL34-mitochondrial ribosomal protein L34 Gene

Reference: PTXBC000071

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L34mt; 39S ribosomal protein L34, mitochondrial; MRP-L34; mitochondrial ribosomal protein L34		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	mitochondrial ribosomal protein L34
<b>GENE ID</b>	64981
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC000071
<b>VECTOR</b>	pDONR233
<b>NCBI SYMBOL</b>	MRPL34
<b>LENGTH</b>	279 bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon )

atggctgtctggctggatccctgttgggcccacgagtaggtcggcagcgttgcgggtggcaggtggctccagccccgggctggctggggttcccagacgcctggggcctccccaccccgagcagggccccgggcaaggctcgcgggaatg  
agtatcagcccagcaacatcaaacgcaagaacaagcacggctgggtccggcgcctgagcacgcggccggcgtgcaggtcatccttcgccgaatgctcaaggcccaagtcgctgagccattga