

MRPL34-mitochondrial ribosomal protein L34 Gene

Reference: PTXBC000071

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L34mt; 39S ribosomal protein L34, mitochondrial; MRP-L34; mitochondrial ribosomal protein L34		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L34
GENE ID	64981
ORIGIN SPECIES	Human
ACCESSIONS	BC000071
VECTOR	pDONR223
NCBI SYMBOL	MRPL34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtcttgctgatccctgttggccccacgtagtgcggcagcgttgcctggcaggtgctccagccccggcctggctgggtccagcgcctgggctccccccccccagcagcggcgggcaaggctcgcgggaatgagatcagccgagcaacatcaa
acgcaagaacaagcagcggctggcgcctgagcagccggcggcgtgcaggtcatccttcgccaatgctcaaggccgaagtcgctgagccattga