

HMGN1-high-mobility group nucleosome binding domain 1 Gene

Reference: PTXBC000075

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG14; non-histone chromosomal protein HMG-14; high mobility group nucleosome-binding domain-containing protein 1; high-mobility group (nonhistone chromosomal) protein 14; high-mobility group nucleosome binding 1; nonhistone chromosomal protein HMG-14; hi		

For research use only.

GENE DESCRIPTION	high-mobility group nucleosome binding domain 1
GENE ID	3150
ORIGIN SPECIES	Human
ACCESSIONS	BC000075
VECTOR	pDONR233
NCBI SYMBOL	HMGN1
LENGTH	303 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagaggaaggtcagctccgccaagggcggccaaggaagagcccaagagagatcggcgggtgtcagctaaacctctgcaaaagtggagcgaagccgaaaaagcagcagcgaaggataatctcagacaaaaaagtca
aacaagaagggaaggagcaaaagggaacaggccgaagtggctaaccaagaactaaagaagacttacctcggaaaacgggaaacgaagactgaggagatccagcctctgatgaagcaggagagaaagaagccaagtctgattaa