

## HMGN1-high-mobility group nucleosome binding domain 1 Gene

Reference: PTXBC000075

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG14; non-histone chromosomal protein HMG-14; high mobility group nucleosome-binding domain-containing protein 1; high-mobility group (nonhistone chromosomal) protein 14; high-mobility group nucleosome binding 1; nonhistone chromosomal protein HMG-14; hi		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	high-mobility group nucleosome binding domain 1
<b>GENE ID</b>	3150
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC000075
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	HMGN1
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atgcccaagaggaaggtcagctccgccgaagcgcgccaaggaagagcccaagagagatcggcgcggtgtcagctaaacctctgcaaaagtggaaagcgaagccgaaaaaggcagcagcgaaggataaatctcagacaaaaagtcaacaaaaggaaaaagggag  
caaaaggaaaacaggccgaagtggctaaccaagaactaaagaagacttacctcggaaaacgggaaacgaagactgaggagagtcagcctgatgaagcagagagaaagaagccaagtctgattaa