

CCND1-cyclin D1 Gene

Reference: PTXBC000076

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCL1; D11S287E; PRAD1; U21B31; G1/S-specific cyclin-D1; B-cell CLL/lymphoma 1; B-cell lymphoma 1 protein; BCL-1 oncogene; PRAD1 oncogene		

For research use only.

GENE DESCRIPTION	cyclin D1
GENE ID	595
ORIGIN SPECIES	Human
ACCESSIONS	BC000076
VECTOR	pDONR223
NCBI SYMBOL	CCND1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggaacaccagctctgtgctgcaagtggaacctccgcccgcgtaccctgatgccaactcctcaacgaccgggtgctgcccggccatgctgaagcggaggagacctgcgcgccctcgggtcctacttcaaatgtgtgcagaaggaggtcctccgtccatgcggaatcgtcgcacactgagaggtctgaggaacagaagtgcgaggaggaggtctcccctgcccctgcccctgaactacctgaccgcttctctgctgctggagccctgaaaaagagccgctgcagctgctggggccactgcatgtctggtccttaagtgaaggaccatccccctgacggcggagaagctgtgcatctacaccgacaactccatccggcccggaggagctgctcaaatggagctgctcctggtgaacaagctcaagtggaaactggccgcaatgaccccgcagatttcaatgaacactcctctccaaaatgccagagggcggaggagaacaacagatcatccgaaaacgcgcagactctggtgcccctgtgcccacagatgtgaagtcaattccaatccgccctccatgggtggcagcggggagcgtggtggccgagtgcaaggcctgaacctgaggagcccaacaactctgtcctactaccgctcacagcttctctccagatgatcaagtggacccgactcctccggcctccagagcagatcgaagccctgctggagtcaagcctgcccaggcccagcagaacatggacccaaggccgaggagaggaagaggaggaggaggtggacctgcttgcacaccaccgactgcccggactggacatcga
```