

RPL8-ribosomal protein L8 Gene

Reference: PTXBC000077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	60S ribosomal protein L8		

For research use only.

GENE DESCRIPTION	ribosomal protein L8
GENE ID	6132
ORIGIN SPECIES	Human
ACCESSIONS	BC000077
VECTOR	pDONR233
NCBI SYMBOL	RPL8
LENGTH	774 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgggccgtgtgatcgtggacagaggaaggcgccgggtctgtgtccgcgcacgtgaagcaccgtaaaggcgtcgcgcctgcgcgccgtggattcgtgagcggcagggctacatcaaggcatcgtcaaggacatccacgacc
cgggcccgcgcccccctcgaagggtgtctccgggatccgtatcggttaagaagcggacggagctgttcattgccgccgagggcattcacacgggccagtttggtaattcggcaagaaggccagctcaacattggcaatgtgtccctgt
gggcaccatgcctgagggtacaatcgtgtgctgctggaggagaagcctggagaccgtggcaagctggccccggccatcagggaaactatgccaccgttatctcccacaacctgagaccaagaagaccgtgtgaagctgccctccggctccaaga
aggttatctctcagccaacagagctgtgtgtgtgtgtgtgtggaggtggccgaattgacaaacccatcttgaaggctggccggcgctaccacaaatataaggcaagaaggaactgctggccacgagtaggggtgtggccatgaatcctgtgga
gcacctttggaggtggcaaccaccagcacatcggcaagccctccaccatccgcagagatgccctgctggccgcaaaagtgggtctcattgctcccgggactggacgtctccgggaaccaagactgtgcaggagaaagagaactag
```