

## COLEC11-collectin sub-family member 11 Gene

Reference: PTXBC000078

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3MC2; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; CLK1; collectin-11; Collectin K1; collectin kidney protein 1; collectin sub-family member 11; collectin subfamily member 11		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	collectin sub-family member 11
<b>GENE ID</b>	78989
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC000078
<b>VECTOR</b>	pDONR233
<b>NCBI SYMBOL</b>	COLEC11
<b>LENGTH</b>	816 bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon )

```
atgagggggaatctgcccctggtggcgttctaatcagcctggccttctgactgctccatctggacatcctcagccggctggcgtgatgacgctgctctgtgcagatcctgctccctggcctcaagggggatgcccggagagagaaggagacaag
gcccccggagggcctggaagagtgcccccacgggagaaaaaggagacatgggggacaaaaggacagaaaaggcagtggtggctgctatggaaaaattggtccattggctctaaaggtgagaaaaggagattccggtgacataggaccccctg
gtcctaattgagaaccaggcctccatgtgagtgacccagctcgcaaggccatcggggagatggacaaccaggctctcagctgaccagcgcagctcaagttcatcaagaatgctgtgccgggtgctgcgagacggagagcaaatctacctg
ctggtgaaggaggagaagcgtctacggcagcccagctgctcctccagggccgccccggcagcctgagcatgcccaaggacgagggctccaatggcctgatggccgatacctggcgcaagccggcctggcccgtcttctatcggcatcaac
gacctggagaggaggcgccttctgtactctgaccactccccatcgggacctcaacaagtggcgcagcgggtgagcccaaatgcctacgacgaggaggactgctggagatggtggcctcggcggtggaacgactggcctgccac
accacatgtacttcatgtgtgagttgacaaggagaacatgtga
```