

COLEC11-collectin sub-family member 11 Gene

Reference: PTXBC000078

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3MC2; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; CLK1; collectin-11; Collectin K1; collectin kidney protein 1; collectin sub-family member 11; collectin subfamily member 11		

For research use only.

GENE DESCRIPTION	collectin sub-family member 11
GENE ID	78989
ORIGIN SPECIES	Human
ACCESSIONS	BC000078
VECTOR	pDONR223
NCBI SYMBOL	COLEC11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggggaatctggccctgtggcgttctaatcagcctggcctctgtcactgtgcatctggacatctcagccggctggcgtgatgacgcctgtctgtcagatctctccttggcctcaaaggggatgctgggagagaagggagacaaaagcgccccggacggcctgaa
gagtcggcccccacgggagaaaaaggagacatgggggacaaaaggacagaaagcagtggtggctgtcatggaaaaatfggtccatggctctaaagtgagaaaaggagattccggtgacataggacccccctgtctctaatggagaaccagcctccatgtgagtcagccagctg
cgcaaggccatcggggagatggacaaccaggctctcagctgaccagcagctcaagttcatcaagaatgctgtcggcgtgtgctgagacggagagcaagatctacctgtgtggaaggagagaagcgtacgggacgccagctgtcctgccaggccgggggacg
ctgagcatgcccaaggacagggctgccaatggcctgatggcgcatacctggcgaagccggcctggcccgtgtctctcatcgacataacgacctggagaaggagggcgcttctgtactctgaccactccccatcgcgacctcaacaagtgccgagcgggtgagccaacaat
gcctacgacgaggaggaactgcgtggagatggtggcctggggcgtggaacgacgtggcctgccacaccatgtactcatgtgtgagttgacaaggagaacatgta