

DEXI-dexamethasone-induced transcript Gene

Reference: PTXBC001083

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Dexi homolog; MYLE; dexamethasone-induced protein; dexamethasone-induced transcript		

For research use only.

GENE DESCRIPTION	dexamethasone-induced transcript
GENE ID	28955
ORIGIN SPECIES	Human
ACCESSIONS	BC001083
VECTOR	pDONR223
NCBI SYMBOL	DEXI
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgccgggtcgcggccacctggacgcactggcccccctgtgccctacgtgccgccgcgtgctgccctctatgtttacgtgggcctgtttcttgcaatgtgctgacctgtactacgccttctcatggagtacatgctcctcaacgtgggcctcgtttctgccccgagg
acatggaccaggcgtcgtggacctggcgtgctcgcgaccccggtcggcctttacgatgctgactcggagctcgcgtctttgatgcgtactggagtag