

ROMO1-reactive oxygen species modulator 1 Gene

Reference: PTXBC008488

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf52; MTGM; MTGMP; bA353C18.2; reactive oxygen species modulator 1; PCM19; ROS modulator 1; epididymis tissue protein Li 175; glyrichin; mitochondrial targeting GXXXG protein; mitochondrial targeting GxxxG motif protein; protein MGR2 homolog		

For research use only.

GENE DESCRIPTION	reactive oxygen species modulator 1
GENE ID	140823
ORIGIN SPECIES	Human
ACCESSIONS	BC008488
VECTOR	pDONR223
NCBI SYMBOL	ROMO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggtggcgtgggtccctacggacagtcagccagccaagctgctcgaccgtgtcaaaatgggctcgtgatgggttgcgcgtgggcagtgccgcccggggcgctcttcggcaccctttctgtctcaggatcggaatgcggggtcgagagctgatggcgccattgggaaacat
gatgcagagtgccgcaccttggcaccatcatggccattgggatgggcatccgatgctaa