

UBC-ubiquitin C Gene

Reference: PTXBC014880

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG20; polyubiquitin-C; ubiquitin C		
For research use only.			

GENE DESCRIPTION	ubiquitin C
GENE ID	7316
ORIGIN SPECIES	Human
ACCESSIONS	BC014880
VECTOR	pDONR223
NCBI SYMBOL	UBC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagatctcgtgaagaccctgactgtaagaccatcactctcaggtggagccgagtgacaccattgagaatgtcaaggcaaatccaagacaaggagccatccctctgatcagcagaggttgatcttctgggaaacagctggaagatggacgacccctgtctgactacaac
atccagaagagtgccaccctgcacctgggtctccgtcttagaggtgggatgcagatctcgtgaagaccctgactgtaagaccatcactctcgaagtgagccgagtgacaccattgagaatgtcaaggcaaatccaagacaaggagccatccctctgaccagcagaggttgatc
tttctgggaaacagctggaagatggacgacccctgtctgactacaacatccagaagagtgccaccctgcacctgggtctccgtcttagaggtgggatgcagatctcgtgaagaccctgactggaagaccatcactctcgaagtgagccgagtgacaccattgagaatgtcaaggcaa
agatccaagacaaggaagcatccctctgaccagcagaggttgatcttctgggaaacagctggaagatggacgacccctgtctgactacaacatccagaagagtgccaccctgcacctgggtctccgtctcagaggtgggatgcagatctcgtgaagaccctgactgtaagacca
tcacctcgaggtggagccagtgacaccatcgagaatgtcaaggcaaatccaagataaggaagcatccctctgatcagcagaggttgatcttctgggaaacagctggaagatggacgacccctgtctgactacaacatccagaagagtgccactctgcacttgggtctcgctt
gaggggggtgtctaa