

MRPS11-mitochondrial ribosomal protein S11 Gene

Reference: PTXBC032378

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCC-2; MRP-S11; S11mt; 28S ribosomal protein S11, mitochondrial; cervical cancer proto-oncogene 2 protein; mitochondrial ribosomal protein S11		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S11
GENE ID	64963
ORIGIN SPECIES	Human
ACCESSIONS	BC032378
VECTOR	pDONR223
NCBI SYMBOL	MRPS11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggctgtgagaacgcgggtcgcgttctcgtcgttggactggccccagacagccggcagggtcgtggccagaacgccggccggaccatctgcacaggcgtcgcagctcgaagacgtcggccaagcagaagtgaacagaacgcggctccagccacacc
aagttcagcatttaccctccattccaggagaggagagctctctgagggtggcaggaaagaaatttgaggagatcccaattgcacacattaaagcatccacaacaacacagatccaggtagtctctgctagttaatgagcccttgccttgcctctgtggcacagaggatttcggaat
gccaagaagggcacaggcatcgcagcacagacagcaggcatagccgcaggcggcagagctaaacaaaaggcgtgatccatccgagttgtgtgaaaggcctggggccaggacgttctgcatgcacggactgatcatggcgccctggaagtgatcctaatcacagaca
acacccaatccacacaacggctgccgccccaggaaggctcggaaagctgtga