



Reference: PTXBC034246

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAS3; CASS3; EFS1; EFS2; HEFS; SIN; embryonal Fyn-associated substrate; Cas scaffolding protein family member 3; signal transduction protein (SH3 containing)		

For research use only.

GENE DESCRIPTION	embryonal Fyn-associated substrate
GENE ID	10278
ORIGIN SPECIES	Human
ACCESSIONS	BC034246
VECTOR	pDONR223
NCBI SYMBOL	EFS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1