

KIAA0895-KIAA0895 Gene

Reference: PTXBC110422

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein KIAA0895		
For research use only.			

GENE DESCRIPTION	KIAA0895
GENE ID	23366
ORIGIN SPECIES	Human
ACCESSIONS	BC110422
VECTOR	pDONR223
NCBI SYMBOL	KIAA0895
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtagccaccggcgtggcgccgctcgggactccacctgcgccccgacccttcggggcgccaggagggcgtcgcgaaaaactgtggagggaagagagaatcggcccttatcttctgctacttattcaaaaagcttcaactggcctgagcaagaactgctaagaagtc
tattctaaatgcagaagattcattgatcattgacaacaaaagacatttcacatttgcctcgggagtgctaaaaagacatttcacaactggaaccagtagttacaatgtcctactacagagcaaggaggaaaaagatcattcacaacacagcttccctccactctccaaaagatgtaga
aaaccagcaaatctcctaactctctgtagcaaaagatcctgcaggatgaaagccctggctcctgtgacaagcagtggtacttgcctgctgagagggcgccctgctgtttttgtcactagttcagtgtaagtcctgtaaaagttacacatgatatctgtttacagggaatggcatagta
ctgccacctaaccacaaaagcaaggtcaagtggtgccatttctccactcttcaagccaaagcctcagctgtctagaagctttgaaaaggagatgactttctgggaagaattttgtatattgactgctataaaaccacacactagagaaagaaaaactgagattcttcaatctgactata
cctacaatcctcagtttgagtatgccaatcctcctctgccaagcgtattagctaaagcatgacccacgcatctgaccgatttcttaagcaggtaaaatgccatttcaatagtggttagttcattttatgtactaaacaattatctgtgacttataagcaaaaacaaagacatttaa