

C19orf30-chromosome 19 open reading frame 30 Gene

Reference: PTXBC117294

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf30; LINC00306; NCRNA00306; PGSF1; MIR7-3 host gene (non-protein coding)		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 30
GENE ID	284424
ORIGIN SPECIES	Human
ACCESSIONS	BC117294
VECTOR	pDONR223
NCBI SYMBOL	C19orf30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgggaatgaggtgtgttcaggttggcgcatggacattttccagaaagggacagagacggcgaagtttgacgtgtctggaagcagagaccagcagggtgactgcttgggagcaccacaaatccggacagcgccctctcgggaggtccgagaagagaaccgcgatctgttca
gcaccggggctcaggacagttccagcgggctcgttttagtctcagaacctggacagctcctccagcttggaaatgcactccctccacctcaaccagagctccccacaactgacctgaccttcttgcagctccatttcatcaagggaacgatccttattgcctcaccatttcccacgt
gaagtctgtattgacattctcatag