

CDH29-cadherin-like 29 Gene

Reference: PTXBC132756

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDH29; PRO34300; cadherin-related family member 4; Cadherin-like protein UNQ9392/PRO34300; VLLR9392; cadherin-like 29; cadherin-like protein 29; cadherin related family member 4		

For research use only.

GENE DESCRIPTION	cadherin-like 29
GENE ID	389118
ORIGIN SPECIES	Human
ACCESSIONS	BC132756
VECTOR	pDONR223
NCBI SYMBOL	CDH29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgctgctcaggctcctcgtgttctctttgtccgggtgtctctgacctctgcagctgccctgtttataaatgtctctgagagccaggccctggcacagctcctcagttttatcctcaactgctcctctacacgccacacccacccctggagttgctcaatgtccagccaccaccac
cttcttcaaccacccagcttggccaggtggcaaggacatgtgggcaagttgacctgagcagctctgctcagttggatgccctgatggtgaaccactacaaggtgcagctgaagttcacatgtggcaacctgtgatggaggctcactctctgtggtgtgcagcgggaccttagcc
atatccagttgtgtcaatttgcagccaggtgagggcaggggcagccggcaggggtggggcaggcatggcctgagcagatcctcctgacctcactctgcccagctggggaaatgattcaggtgccagagacagtcacacctggggctcggctgtacactctgctcctccag
gcctagaactccacggagcaggttaagccagaacatgggacgggtgggtga