

**FAM101A-family with sequence similarity 101, member A Gene**

Reference: PTXBC141805

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM101A; filamin-interacting protein FAM101A; FAM101A; CFM2; family with sequence similarity 101, member A; refilinA; regulator of filamin protein A; refilin A		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	family with sequence similarity 101, member A
<b>GENE ID</b>	144347
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC141805
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	FAM101A
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

```
atgaggccccggatgctgccagtgttcttggggagagcatcaaggtgaaccggaaaccacgcatgagatccgctgcaactctgaggtcaagtaagcctcggagaagcattccaggacaaggtcttctatgcgctgtaccaccgtcaaggcctacagcagaccatcgtggcagc
accctcaactgcagctggcgcaactaccgagccagctgacctggagccacgcccgcgcgacctgcctccgcagcaccaccatcatctcccgaagcatgccaggagcaacttccggaccaccctgcaactgcagcctgggcccggcccagccgctggttaccgccagcctgcagct
gcagcttggcaggaccctgccccagcctctgggacctgccacgctctga
```