

CNN2-calponin 2 Gene

Reference: PTXBC141818

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	calponin-2; calponin H2, smooth muscle; neutral calponin; calponin 2		

For research use only.

GENE DESCRIPTION	calponin 2
GENE ID	1265
ORIGIN SPECIES	Human
ACCESSIONS	BC141818
VECTOR	pDONR233
NCBI SYMBOL	CNN2
LENGTH	930 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgagctccacgcagttcaacaaggccctctgtacggctgtcggccgaggtcaagaaccgctctgtcacaatatgacccccagaaggagcagagctccgcacctggatcagggaactcaccggcctccatggccccgaactccaagaaggcctgaaggatggaactatc
ttatgcacactcatgaacaagctacagccgggctccgtccccaagatcaaccgctccatgcagaactggcaccagctagaaaacctgtccaactcacaaggccatggtcagctacggcatgaaacctgtggacctgtcgaaggccaacgacctgttgagatgggaacatgacgcag
gtgcaggtgtctctctcgcctggctgggaaaggccaaactaaagggtcgaagcggggtggacattggcgtcaagtactcgaagaagcaggaagcggaaattcgaacgatgccaccatgaaggctggccaagtgcgtcaccggctgcaatggcaccacaacaaatgcgccagcca
gtcgggcatgactgcctacggcacgagaaggcatctatgaccccaagaaccatctctgcccccatggaccactcaccatcagcctccagatgggcacgaacaagtgccagccaggtgggcatgacggctcccgggaccggcgccaca tctatgatacgaactgggaac
cgacaaagtgtgacaactctccatgtccctgcagatgggtacacgagggcgccaaccagagcggccaggtcttcggcctggccggcagatatatgaccccaagtactgcccgaaggcacagtgccgatggggtccctcgggcaccggcgactgccccgacccgggggag
gtccctgaatatcccccttactaccagaggagggcgtactga
```