

FCGR2B-Fc fragment of IgG, low affinity IIb, receptor (CD32) Gene

Reference: PTXBC141826

| | | | |
|----------------|---|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | CD32; CD32B; FCG2; FCGR2; IGFR2; low affinity immunoglobulin gamma Fc region receptor II-b; CDw32; Fc fragment of IgG, low affinity II, receptor for (CD32); Fc fragment of IgG, low affinity IIb, receptor (CD32); Fc fragment of IgG, low affinity IIb, recep | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | Fc fragment of IgG, low affinity IIb, receptor (CD32) |
| GENE ID | 2213 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC141826 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | FCGR2B |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atgggaatctgtcattctacctgtccttgccactgagagtactggctgactgcaagtccccagccttggggtcatatgcttctgtggacagctgtgctattctggctctgttctgggacacctgcagctccccaaaggctgtgctgaaactcgagccccagtggatcaactgtctccaggaggactctgtgactctgacatccgggggactcacagcctgagagcgactccattcagtggttcacaatgggaatctcattcccaccacacgcagcccagctacaggttcaaggccaacaacatgacagcggggagtacacgtgccagactggccagaccagctcagcgacctgtgcatctgactgtccttctgagtggtgctcctcagaccctcactggagttccagggaggagaacctcgtgctgaggtgccacagctggaaggacaagcctctggtcaaggctacattctccagaatggaaaatccaagaaatctccgttggatcccaactctccatcccacaagcaaacacagctcacagtggtgattaccactgcacaggaacataggtacacgctgtaactatccaagcctgtgaccatcactgtccaagctccagctctcaccgatgggatcattgtgctgtgtgctactgggactgtgtagcggcattgttctgctgtagtggtgatctactgcaggaaaaagcggatttcagcaatcccactaatcctgatgaggtgacaagttgggctgagacaacaatcacctattcaactctcatgcccggatgctctggaagacctgatgaccagaacctatttag