

**FCGR2B-Fc fragment of IgG, low affinity IIb, receptor (CD32) Gene**

Reference: PTXBC141826

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD32; CD32B; FCG2; FCGR2; IGFR2; low affinity immunoglobulin gamma Fc region receptor II-b; CDw32; Fc fragment of IgG, low affinity II, receptor for (CD32); Fc fragment of IgG, low affinity IIb, receptor (CD32); Fc fragment of IgG, low affinity IIb, recep		

For research use only.

---

<b>GENE DESCRIPTION</b>	Fc fragment of IgG, low affinity IIb, receptor (CD32)
<b>GENE ID</b>	2213
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC141826
<b>VECTOR</b>	pDONR233
<b>NCBI SYMBOL</b>	FCGR2B
<b>LENGTH</b>	876 bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

```
atgggaatcctgtcattcttacctgtccttgccactgagagfgactgggctgactgcaagfccccccagccttggggtcatatgcttctgtggacagctgtgctattcctggctcctgttgcctgggacacctgcagctcccccaaaggctgtgctgaaactcgaagccccagtgatcaacgtgct
ccaggaggactctgtgactctgacatgccgggggactcacagccctgagagcgaactccattcagtggtccacaatgggaatctcattccaccacacgcagcccaagctacaggttcaaggccaacaacaatgacagcgggggagtacacgtgccagactggccaagaccagcctcagc
gacctgtgcatctgactgtgctttctgagtgctgtgctccagaccctcaccctggagttccaggaggagaaaccatctgtctgaggtgccacagctggaaggacaagcctctgtcaaggtcacattctccagaatggaaaatccaagaaatttcccgttggatcccaacttctca
tcccacaagcaaacacagtcacagtggtgattaccactgcacagaaacataggctacacgctgtactatccaagcctgtgacacatcttccaagctcccagcttccaccgatgggacattgtggctgtgctcactgggactgctgtagcgccattgttgcgtgtatggcctt
gatctactcagaaaaagcggatttcagccaatcccactaatctgataggctgacaaaagttgggctgagaacacatcactattcctctcatgcacccggatgctctggaagacctgatgaccagaaccgtatttag
```