

RPS15-ribosomal protein S15 Gene

Reference: PTXBC141832

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RIG; S15; 40S ribosomal protein S15; homolog of rat insulinoma; insulinoma protein; ribosomal protein S15		

For research use only.

GENE DESCRIPTION	ribosomal protein S15
GENE ID	6209
ORIGIN SPECIES	Human
ACCESSIONS	BC141832
VECTOR	pDONR233
NCBI SYMBOL	RPS15
LENGTH	438 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggcagaagtagcagaaagaagcggacctccgcaagttcacctaccggcgtggacctgaccagctgctggacatgctctacgagcagctgatgcaagctgtacagtgccgccagcggcggcggctgaaccggggcctcggcggaaagcagcactccctgctgaagc
gctgcgcaaggccaagaaggcggcggccatgggaagccggaagtggfgaagacgcacctggggacatgatcatctaccgagatggfggcagcatggggcgtctacaacggcaagacctcaaccaggtggagatcaagcccagatgctggccactacctgg
gcgagttctccatcacctacaagcccgtaaagcatggcggccggcatcggggccaccactctcccgttcatccctcaagtaa
```