

GOLGA2L1-golgi autoantigen, golgin subfamily a, 2-like 1 Gene

Reference: PTXBC141930

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	golgi autoantigen, golgin subfamily a, 2-like 1
GENE ID	55592
ORIGIN SPECIES	Human
ACCESSIONS	BC141930
VECTOR	pDONR233
NCBI SYMBOL	GOLGA2L1
LENGTH	435 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacatgagggaggaggaggaggtgcctcagcccatgccaagcatcccggaggatctagagagccagaaggccatggtggcattttcaactcagctgtagctagtccgaggaggagcagcagcggctatgtggcagctgaaggagtgcactgccagcgcctggctcatctgtggcctcggcccaagaaggaaacctgagcagcagccccagccccagaactgggggatccatgtgtgggagaccaccagggccctgcagggggccatggagaagctgtgggagagtacatgcactgtaccagagccagggcagctcgggaaggaggaggtgcatcagcagctgccccagacaaggaggaggtgaaagtgaagctgctggagctggcgtggccttggagactgcaacaagtggcatag