

LOC646938-similar to TBC1 domain family, member 2B Gene

Reference: PTXBC141940

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	similar to TBC1 domain family, member 2B
GENE ID	646938
ORIGIN SPECIES	Human
ACCESSIONS	BC141940
VECTOR	pDONR233
NCBI SYMBOL	LOC646938
LENGTH	258 bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacaggcaggaaagtcctcgggcatctctgttcaacacagctttctgggctggaagctgacagatctctttgggacctgaaaccttccctacgccagatccggaaccaacgcgcctaccacttgagaaaagtcggctggagctgaccgagctggaggccatccg
tgaggacttctgcatgagcgggacaccgccctgacaagggtgagctgtcagtgacgaggaggagatacctga