

E. coli ArsR Nter GST and FLAG Recombinant Protein

PX-P1070-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	E. coli
FRAGMENT	Partial (Δ Cter)
PROTEIN DELIVERED WITH TAG	No
MOLECULAR WEIGHT WITH TAG IF ANY	14,55 kDa
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	0,8
PROTEIN ACCESSION	ZP_03049476.1
FORM	liquid
BUFFER	TrisHC 50mMl, 10mM reduced glutathion pH8 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

GPLGSDYKDDDDKGSFLLPIQLFKILADETRLGIVLLLSELGELCVCDLCTALDQSQPKISRHLALLRESGLLLDRKQGGK
WVHYRLSPHIPAWAAKIIDEAWRCEQEKVQAIVRNLARQNCSGDSKNICS

More info

SPECIFICATION

GENE ID
SWISSPROTID
NCBI GENE ALIASES
UNIPROT ID
UNIPROT LINK

SYNONYMS

ArsR, arsenical resistance operon repressor