

## Staphylococcus GlmCat Recombinant Protein

PX-P2072-10

### OVERVIEW

---

<b>SIZE</b>	10 ug
<b>ORIGIN SPECIES</b>	Staphylococcus
<b>FRAGMENT</b>	
<b>PROTEIN DELIVERED WITH TAG</b>	Yes
<b>MOLECULAR WEIGHT WITH TAG IF ANY</b>	38.20 kDa
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	E.coli
<b>PURITY</b>	90%
<b>PROTEIN ACCESSION</b>	KEK36855.1
<b>FORM</b>	Frozen
<b>BUFFER</b>	PBS, pH 7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

### PROTEIN SEQUENCE

---

MAHNHRHKHKL VNAKDLTAPTAVKPTTSAAKDYNITYVIKNGNGYVTPNSDTAKYSLKAFNEQPFVAVKEQVINGQTWYYGKLSNGKLA WIKSTDLAKEL  
IKYNQTGMTLNQVAQIQAGLQYKPQVQRVPGKWTGANFNDVKHAMDTKRLAQDPALKYQFLRLDQPQNISIDKINQFLKGGVLENQGA AFNKA AQMYGINE  
VYLISHALLETGNGTSQLAKGADV VNNKVV TNSNTKYHNVFGIAAYDNDPL

### More info

#### SPECIFICATION

**GENE ID**  
5330556  
**SWISSPROTID**  
A6QFR2  
**NCBI GENE ALIASES**  
**UNIPROT ID**  
**UNIPROT LINK**

#### SYNONYMS

GlmCat, Bifunctional  
autolysin