

Human CD9 Recombinant Protein

PX-P3040-10

OVERVIEW

SIZE	10 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG	10.5 kDa
IF ANY	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	E.coli
PURITY	90%
PROTEIN ACCESSION	
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

PROTEIN SEQUENCE

SHKDEVIKEVQEFYKDTYNKLKTKDEPQRETLKAIHYALNCCGLAGGVEQFISDICPKKDVLETFVTKSCPDAIKEVFDNKFHI

More info

SPECIFICATION

GENE ID

SWISSPROTID

NCBI GENE ALIASES

UNIPROT ID

P21926

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P21926>

SYNONYMS

CD9; TSPAN29; 5H9;
BA2; BTCC-1; DRAP-27;
GIG2; MIC3; P24;
Tetraspanin 29; Cell
growth-inhibiting gene 2
protein; Leukocyte antigen
MIC3; Motility Related
Protein