

Human CD278/ICOS recombinant protein

PX-P4035-100

OVERVIEW

SIZE	100 ug
ORIGIN SPECIES	Human
FRAGMENT	
PROTEIN DELIVERED WITH TAG	Yes
MOLECULAR WEIGHT WITH TAG IF ANY	42.34kDa
DELIVERY CONDITION	any condition

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	mammalian
PURITY	0,69999999999999996
PROTEIN ACCESSION	
FORM	Lyophilized
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

PROTEIN SEQUENCE

MKSGLWYFFLCLRIKVLGTGEINGSANYEMFIFHNGGVQILCKYPDIVQQFKMQLLKGGQILC
 DLTKTKGSGNTVSIKSLKFCHSQLSNNSVSFFLYNLDSHANYYFCNLSIFDPPPFKVTLTGGY
 LHIYESQLCCQLKSRDDDDKDKTHTCPPCPAPPELLGGPSVFLFPPKPKDQLMISRTPEVTCVTV
 DVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKV
 NKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPE
 NNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSVLHEALHNHYTQKSLSLSPGK

More info

UNIPROT ID

Q9Y6W8

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q9Y6W8>

SYNONYMS

CD278, AILIM, Inducible
 T-Cell Costimulator,
 Activation-Inducible
 Lymphocyte
 Immunomediatory
 Molecule

For research use only.